


Extrait de la notice ROCHE ACCU CHECK INSIGHT

7.2 Date, heure et format de l'heure

Menu principal > Réglages > Date et heure


1



The screenshot shows the main menu of the application. At the top, the status bar displays '09:29' and '29 mar 2014'. The menu consists of several icons: 'Lecteur', 'Bolus', 'Communication', 'Mes données', 'Mode avion', and 'Réglages'. The 'Réglages' icon, which shows a gear, is highlighted with a blue background. At the bottom, there are two buttons: 'Etat' and 'Mode utilisateur'.

Dans le menu principal, sélectionnez **Réglages**.

2



The screenshot shows the 'Réglages système' (System Settings) menu. The title is 'Réglages système'. The menu items are: 'Rappels', 'Date et heure', 'Conseil bolus', 'Plages horaires', 'Réglages mode', and 'Unité de glucides'. The 'Date et heure' item is selected and highlighted in blue. At the bottom, there is a back arrow and the word 'Retour'.

Sélectionnez **Date et heure**.

3



The screenshot shows the 'Date et heure' (Date and Time) settings screen. At the top, it says 'Date et heure'. Below that, there are three fields for the date: '29' (Jour), 'mar' (Mois), and '14' (Année). Below the date fields, there are two fields for the time: '09' (h) and '29' (min). Below the time fields, there are two radio buttons: '12 h' and '24 h'. Below the radio buttons, there is a large black button with a white '+' sign and a white '-' sign. At the bottom, there are two buttons: 'Annuler' and 'Enregistrer'.

Date et heure : sélectionnez et définissez les champs appropriés. **Format de l'heure** : choisissez **12 h** ou **24 h**. Sélectionnez **Enregistrer**.